

# Watercare Laboratory Services

## Bottle Guide (Solids)

# Watercare

## Laboratory Services

### General Chemistry

SI1. 250mL Plastic Pottle

SI2. 500mL Plastic Pottle

**Used for:** General chemistry testing such pH and total nitrogen as well as metals testing

**Sampling info:** Please use clean, non-metallic equipment to transfer solid material into the pottle. Ensure the lid is correctly seated on the pottle and screw down firmly. Note, larger volume of samples will be needed for more comprehensive testing.

**Transport info:** Ideally, samples should be transported immediately to the lab for analysis.

**Label:**

Watercare Laboratory Services		Sample
Date Sampled		
Sampled By		
Sample #		
Sample Location		
Client Preparation		
BOTTLES SUPPLIED BY WATERCARE LAB SERVICES E-mail: clientsupport@water.co.nz Ph: (09) 539 7600 52 Aintree Ave, Airport Oaks, Auckland		



### Organic Chemistry

SO1. 250mL Glass

SO2. 500mL Glass

SO3. 1L Glass

SI1. 250mL Plastic Pottle (for PFAS analysis)

**Used for:** Organic chemistry testing such TPH, SVOC, VOC etc.

**Sampling info:** Please use clean, metallic equipment to transfer solid material into the pottle. For PFAS analysis, plastic containers must be used (see SI1). Ensure the lid is correctly seated on the bottle and screw down firmly. Note, larger volume of samples will be needed for more comprehensive testing.

**Transport info:** Ideally, samples should be transported immediately to the lab for analysis.

**Label:**

Watercare Laboratory Services		Sample
Date Sampled		
Sampled By		
Sample #		
Sample Location		
Client Preparation		
BOTTLES SUPPLIED BY WATERCARE LAB SERVICES E-mail: clientsupport@water.co.nz Ph: (09) 539 7600 52 Aintree Ave, Airport Oaks, Auckland		



## Microbiology

### SM1. 400mL Sterile Plastic

### SM2. 1L Sterile Plastic

**Used for:** Microbiology testing such as E. coli, Campylobacter, Viruses etc.

**Sampling info:** Please use clean, sterile equipment to transfer solid material into the pottle. Ensure the lid is correctly seated on the pottle and screw down firmly. Note, larger volume of samples will be needed for more comprehensive testing.

**Transport info:** Samples must be transported immediately and must arrive at the lab within 24 hours of sampling. Keep samples below 10°C.

**Label:**

Watercare Laboratory Services	Sample
Date Sampled .....	
Sampled By .....	
Sample # .....	
Sample Location .....	
Client Preparation .....	
BOTTLES SUPPLIED BY WATERCARE LAB SERVICES E-mail: <a href="mailto:clientsupport@water.co.nz">clientsupport@water.co.nz</a> Ph: (09) 539 7600 52 Aintree Ave, Airport Oaks, Auckland	

